

Accessing the Genomic Dark Matter: Using Next-Generation Mapping to Identify Large-Scale Repeat Structures in Maize

Jinghua Shi¹, Xiang Zhou¹, R. Kelly Dawe², Ahmed Naguib¹, Zhanyang Zhu¹, Zeljko Dzakula¹, Weiping Wang¹, Saki Chan¹, Alex Hastie¹, Han Cao^{#1}

¹ BioNANO Genomics, Inc., San Diego, CA

² Department of Plant Biology, Athens, University of Georgia, GA


Abstract

Complex genome assembly has improved dramatically with advent of new long-range sequencing and assembly technologies. Among these methods, the next-generation mapping (NGM) with BioNANO Irys® system offers distinct advantages by generating optical maps using NanoChannel arrays technology, particularly when combined with long range sequencing approaches such as PacBio.





The recently released maize B73 reference genome was assembled using a hybrid scaffolding approach that combined PacBio sequence and BioNANO maps. While this is an excellent assembly with long contig N50, there are areas where the BioNANO map is continuous but matching sequence data could not be identified. We analyzed BioNANO maps corresponding to the largest of these sequence gaps (>100 kbp) and found that many are associated with either large duplicated regions or short tandem repeat arrays. These include the 45S and 5S rDNA arrays and several centromeric repeat arrays. In addition, we show that the P1 locus, which has 11 nearly-perfect copies in the B73 reference genome, has 10 copies in another maize inbred W22, demonstrating the utility of using NGM to assess copy number variations (CNV).

The data also reveal multiple previously unknown large-scale CNVs, many of which likely include genes. These results point to limitation of current sequencing methodologies and highlight the value of using NGM generated ultra-long 150 kbp -2 Mbp molecules for the accurate assembly of the maize genome, a process that is both cost effective and fast.

Methods



Results



Conclusions

The BioNANO Genomics Irys System provides a platform for mapping and visualizing extremely long molecules, which makes it possible to analyze long tandem repeat arrays comprehensively, accurately and without bias due to selection or amplification. Accompanied by our repeat detection tool (X. Zhou et al. PAG XXIV poster), repetitive units of the complex genomes can be profiled, analyzed and quantified in an unprecedented fashion. Our results show that tandem repeat arrays can be reliably detected as regularly spaced labeling patterns (as is the case for the P1 locus), or as long no label regions (as is the case for centromere 10) provided there are diagnostic patterns on either side of the array.

See also Posters: P0957, P0961, P0033, and P0958.

References

- Wolfgruber, T. K., M. M. Nakashima, K. L. Schneider, A. Sharma, Z. Xie et al., 2016. Frontiers in Plant Science 7.
- Springer, N. M., K. Ying, Y. Fu, T. Ji, C. T. Yeh et al., 2009. PLoS Genet 5: e1000734.
- Coccilone, S. M., D. Nettleton, M. E. Snook and T. Peterson, 2005. Plant Biotechnol J 3: 225-235.
- Jiao, Y., P. Peluso, J. Shi, T. Liang, et al. in revision.
- X. Zhou, S. Chan, J. Shi, et al. PAG XXIV poster.