

Jen Hauenstein, Andy Pang, Alka Chaubey, Alex Hastie

Bionano Genomics, Inc., San Diego, CA 92121, USA  
Correspondence: ahastie@bionanogenomics.com


Abstract

Guidelines from WHO, NCCN and others, for the genetic analysis of hematological malignancies included structural variation analysis. Traditionally, this has relied on a combination of three cytogenetic technologies for structural variation analysis: karyotyping, FISH, and microarray, are used to detect copy number variants, translocations, inversions. Next generation sequencing is also applied for mutation analysis but has not been successful for structural variation analysis. These traditional methods have many very manual aspects and require extensive expertise. Optical genome mapping (OGM) consolidates assays into a single laboratory assay in which the output provides the visualization of structural and copy number variants at one time.

OGM is able to comprehensively detect structural variations genome wide down to 5% variant allele fraction for CNVs, inversions, and translocations from blood and bone marrow aspirates making it an attractive choice for hematologic

malignancy genomic analysis. Preanalytical and analytical steps require approximately 4-5 days from sample to processed data with structural variation calls. Dynamic filtering in the user interface can be configured to remove most polymorphic variants and prioritize relevant variants. In addition, the OGM graphical user interface software, Bionano Access 1.7, allows for the user to assign classification/relevance to the variants for each case. For example, an ALL sample with t(9;22), deletion of *CDKN2A*, and whole chromosome gains of 4,6, and 10 can be easily visualized with the Circos plot and, then, can be further examined and annotated as needed. A second analyst can repeat the process blind to the first analysis and a supervisor can adjudicate the classifications. A variety of cases with hallmark abnormalities from various leukemias will be presented with the filtering and prioritization workflow used to derive them. This comprehensive technology allows for a quicker, more reliable output than traditional cytogenetic approaches.

Methods




(1) Long molecules of DNA are labeled with Bionano reagents by (2) incorporation of fluorophores at a specific sequence motif throughout the genome. (3) The labeled genomic DNA is then linearized in the Saphyr chip using NanoChannel arrays (4) Single molecules are imaged by the Saphyr instrument and then digitized. (5) Molecules are uniquely identifiable by distinct distribution of sequence motif labels (6) and then assembled by pairwise alignment into *de novo* genome maps.


## Virtual Panels

All classes of structural variations can be detected by OGM in a single assay. Above, deletions, amplifications, inversions, and translocations are shown to represent all variants recommended by the National Comprehensive Cancer Network (NCCN).


### Acute Myelogenous Leukemia



### Chronic Myelogenous Leukemia



### Acute Lymphocytic Leukemia



Not shown: CEP4, 10, 17, 9p21/9q21

## Publications

Reference	Cohort size	Clinical referral	Samples type	Variant types	Number of abnormalities included	Concordance with cytogenetics results	OGM additional findings
Radbound University Nevling et al., 2020	48	Myeloid and lymphoid neoplasms (AML, MDS, CML, CLL, ALL, MLL, MPN, T-PLL, LYM)	Peripheral blood or bone marrow	Copy number gains/losses, inversions, translocations, aneuploidies, LOH	112	100%	18 potential gene fusions absent from COSMIC database 26 insertions/deletions overlapping with well-established cancer genes
Cancer Genomics Consortium Ley et al., 2020	100	Acute myeloid leukemia	Peripheral blood or bone marrow	Copy number gains/losses, inversions, translocations, aneuploidies	NA	100%	3 translocations, 1 inversion, 2 deletions and 1 derivative chromosome
CHU Amlens Lestragant et al., 2021	10	B and T Acute lymphoblastic leukemia	Peripheral blood or bone marrow	Copy number gains/losses, inversions, translocations, aneuploidies, LOH	78	78 out of 78*	4 fusions, 0 deletions, 2 gains, 1 duplication, 3 complex chromosomal rearrangements
M.D. Anderson Tang et al., 2020	12	Myelodysplastic Syndromes	Fresh/ frozen bone marrow aspirate	Copy number gains/losses, inversion, aneuploidies, translocations/ derivative chromosomes and 1 isodicentric chromosome	28	26 out of 28** 93% with conventional karyotyping, 100% with chromosomal microarray	8 cryptic aberrations (in 33% patients), 4 deletions, 1 duplication and 1 translocation
Johns Hopkins University Siretti et al., 2021	5	Leukemia/Lymphoma and Solid Tumors	Peripheral blood, bone marrow, or fresh kidney tissue	Copy number gains/losses, inversions, translocations, aneuploidies	30	100% CMA >10% VAF ***	71 additional calls (7.7% involving cancer genes)
University Hospital Olomouc Kriegova et al., 2021	11	Multiple myeloma	Bone marrow	Copy number gains/losses, inversions, translocations, aneuploidies	NA	98%	

## Conclusions

- ✓ OGM allows for the detection of all classes of structural variations occurring in the human genome
- ✓ Removes the need for tiered testing as part of a clinical work-up
- ✓ OGM can pick up all actionable SVs as per NCCN/WHO guidelines
- ✓ Bionano's Access software solution - ease of use
- ✓ Graphical user interface allows filtering and manual evaluation on a per SV basis including access to raw data support
- ✓ Fast TAT (sample to answer in 4 days)
- ✓ Multiple studies have demonstrated 100% concordance with cytogenetic methods (karyotype, microarray, FISH) for constitutional and hematologic malignancy
- ✓ Multiple studies showing increase in diagnostic yield against the standard of care in pediatric and heme